

Verbiage for ID For Life Members' Website

- **OralID®** is manufactured by Forward Science, a company that was founded with the vision to create a simpler, efficient, and more affordable oral cancer screening device. **OralID®** uses the same proven fluorescence technology as other oral exam devices, with no up-front costs and lower overall lifetime costs. In addition, unlike any other device on the market, **OralID®** has no per-patient cost during use.
- **OralID®** uses a proven, optically based technology called, "fluorescence technology." **OralID®** fluorescence technology uses a blue light (435–460 nm) that provides the clinician an aid in visualization of oral mucosal abnormalities, such as oral cancer and pre-cancer. The other medical applications of fluorescence technology include highlighting abnormalities of the GI tract, cervix, lung, and skin.
- **OralID®** is an adjunctive screening device that uses fluorescence technology to aid clinicians in visualization of oral abnormalities. It is easy to use and non-threatening to patients. Dentists and hygienists can use **OralID®** to increase the efficacy of oral exams without incurring additional per-patient costs.